


THE VIAFLO PLASTIC CONTAINER HANDLING INSTRUCTIONS




- The **VIAFLO Container product is fragile and may be easily damaged if dropped or improperly handled**
- The VIAFLO Container product and associated cartons should be handled with care to avoid damage and the potential for leakage
- The VIAFLO Containers should be **distributed in their original Baxter shipper cartons** during the Customer supply-chain



- When opening the VIAFLO Container shipper carton do not use sharp instruments that may cut the VIAFLO Containers inside the shipper carton



- Do not stack the individual VIAFLO Container shipping cartons more than eight (8) high



- Do not double stack the VIAFLO Container shipping pallets

Product Code	NDC Number	Product Code Description	Pack Factor
UE1322D	0338-9543-02	0.9% Sodium Chloride Injection USP in VIAFLO plastic container, 250 mL	30
FE1323D*	0338-9542-05	0.9% Sodium Chloride Intravenous Infusion BP in VIAFLO Container, 500mL DAC	20
UE1323D	0338-9543-04	0.9% Sodium Chloride Injection USP in VIAFLO plastic container, 500 mL	20
UE1324D	0338-9543-06	0.9% Sodium Chloride Injection USP in VIAFLO pastic container, 1000 mL	10

* The FDA has allowed the importation and distribution of FE product codes listed above.

www.baxter.com

Baxter Healthcare Corporation
One Baxter Parkway
Deerfield, IL 60015

Baxter and Viaflo are trademarks of Baxter International Inc.

